


Taicron Pale Shade Tricombination

Taicron Yellow E-3GT 200%

Taicron Yellow E - 3GT 200%

Taicron Red E- FBLN-T 200%



Taicron Blue E- RNT 200%

Taigal AF-60 0.3 g/l



Dayacid HP PH 4.0

130°C X 30min



Taicron Red E - FBLNT 200%

Taicron Blue E - RNT 200%

Taicron Pale Shade Tricombination

Taicron Yellow E-3GT 200%

Taicron Yellow E - 3GT 200%

Taicron Red E- FBLN-T 200%

0.3	
0	0

Taicron Blue E- RNT 200%

Taigal AF-60 0.3 g/l

0.25		0.25	
0.05	0	0	0.05

Dayacid HP PH 4.0

130°C X 30min

0.2		0.2		0.2	
0.1	0	0.05	0.05	0	0.1

0.15		0.15		0.15		0.15	
0.15	0	0.1	0.05	0.05	0.1	0	0.15

0.1		0.1		0.1		0.1		0.1	
0.2	0	0.15	0.05	0.1	0.1	0.05	0.15	0	0.2

0.05		0.05		0.05		0.05		0.05		0.05	
0.25	0	0.2	0.05	0.15	0.1	0.1	0.15	0.05	0.25	0	0.25

0		0		0		0		0		0		0	
0.3	0	0.25	0.05	0.2	0.1	0.15	0.15	0.1	0.2	0.05	0.25	0	0.3

Taicron Red E - FBLNT 200%

Taicron Blue E - RNT 200%